

Touchless Thermometer Quick Start Instructions

1. Install Batteries
2. Squeeze Trigger until thermometer powers on
3. Quickly hold large side button for a few seconds until **F1** is displayed
4. Press large side button twice to toggle to **F3**
5. Press either of the small buttons to toggle temp **from C to F**
6. Press large side button to exit setup

Using the same method as above, you can also set:

F1: Temperature method. Object (Image of house) or Human (Smile Icon)

F2: Sound. Toggle on/off

F3: Temperature display. Toggle from Celsius to Fahrenheit

TO USE:

Hold close to forehead.

Squeeze trigger and release.

Temperature will be displayed.

Keep in mind, touchless thermometers are not as accurate as in-ear thermometers.

Touchless Thermometer Quick Start Instructions

1. Install Batteries
2. Squeeze Trigger until thermometer powers on
3. Quickly hold large side button for a few seconds until **F1** is displayed
4. Press large side button twice to toggle to **F3**
5. Press either of the small buttons to toggle temp **from C to F**
6. Press large side button to exit setup

Using the same method as above, you can also set:

F1: Temperature method. Object (Image of house) or Human (Smile Icon)

F2: Sound. Toggle on/off

F3: Temperature display. Toggle from Celsius to Fahrenheit

TO USE:

Hold close to forehead.

Squeeze trigger and release.

Temperature will be displayed.

Keep in mind, touchless thermometers are not as accurate as in-ear thermometers.